



HEKA
a division of Harvard Bioscience, Inc.



UPGRADE your EPC10 Patch Clamp Amplifier to an

EPC10 USB



BENEFITS

- Your EPC10 USB amplifier will be compatible by the latest software and the next version. The next version will not support ITC-1600 boards. Please get in touch with us for details.
- The EPC10 USB does not require any driver installation
- The EPC10 USB uses a USB2.0 port, no PCI-card required
- The EPC10 USB can be used with laptop as well as desktop computers
- One year warranty on the new interface

The upgrade will involve replacing the ITC-1600 interface with the current LIH 8+8 interface. Further information about the capabilities of the LIH 8+8 interface can be found at:

www.heka.com/products/products_main.html#acq_lih88

The amount of channels, dynamic range and sample rate remains unchanged.

To qualify for the upgrade your amplifier must be version "J" or later. To verify the version, start PATCHMASTER and go to the EPC10 amplifier menu as demonstrated in the screen shot.

Contact HEKA

HEKA Instruments Inc.
Toll free: +1 800 272 2775
Phone: +1 508 893 8999
Fax: +1 508 429 5732

HEKA Elektronik
Dr. Schulze GmbH
Phone: +49 (0)6325 9553-0
Fax: +49 (0)6325 9553-50