



Pre-amplifiers of the EPC 10 patch clamp amplifiers are available with two different mounting plates for mechanical connection of the pre-amplifier to e.g. various micromanipulator system:

- **Standard Mounting Plate (GZ0012).** The headstage is mounted on a 38 mm x 90 mm x 4 mm plate which has 4 wholes with 3 mm diameter. Since the plate is wider than the headstage there is room for custom mounting wholes on both sides of the headstage.



- **Dovetail Mounting Plate (GZ0013).** The dovetail adapter can be used to connect e.g. to Sutter MP 285 or HEKA MIM 4 micromanipulator. The dimensions of the dovetail match the headstage housing.

