



PSA 12 Tone Generator

Membrane Resistance Monitor for PATCHMASTER



If you believe your job is to
make new discoveries...

...discover PSA 12:

- Sound support for any amplifier connected to PATCHMASTER
- Plug and Play connection via USB 1.1 / 2.0

HEKA

HEKA provides the finest instruments today
to achieve the needed progress of tomorrow...



PSA 12 Tone Generator

Membrane Resistance Monitor for PATCHMASTER

Introduction

The PSA 12 is a USB tone generator which allows amplifiers without built-in sound support to take full advantage of the membrane resistance monitor feature of PATCHMASTER.

With the PSA 12 the PATCHMASTER software has the capability to play a sound that is proportional to R-membrane. Which is typically used for monitoring seal-resistance. The PSA 12 has a frequency range between 100 Hz and 4 kHz. Both the tone characteristics and volume are controlled by the PATCHMASTER software. A standard 3.5 mm audio connector provides connectivity to an optional amplified speaker or to the supplied headset.

Features

- Frequency range between 100 Hz and 4 kHz
- Connectivity to external speakers via standard audio connector
- Software control of volume and tone by PATCHMASTER
- Fully compliant with USB 1.1 / 2.0 specification
- Self powered device
- Includes USB cable and headphones
- Supported by MacOS X (OS 10.39 or higher) and Windows XP / Vista *

* Compatible with Windows XP, Vista and MacOS X (10.39 or higher) without an additional driver installation. Requires PATCHMASTER version 2.15 or higher, with PSA 12 support enabled in the configuration settings.



HEKA Elektronik
Dr. Schulze GmbH
Wiesenstraße 71
D-67466 Lambrecht/Pfalz
Germany

Phone +49 (0) 63 25 / 95 53-0
Fax +49 (0) 63 25 / 95 53-50
Web Site <http://www.heka.com>
Email sales@heka.com
support@heka.com

HEKA Electronics Incorporated
47 Keddy Bridge Road
R.R. #2
Mahone Bay, NS B0J 2E0
Canada

Phone +1 902 624 0606
Fax +1 902 624 0310
Web Site <http://www.heka.com>
Email nasales@heka.com
support@heka.com

HEKA Instruments Inc.
2128 Bellmore Avenue
Bellmore, New York 11710-5606
USA

Phone +1 516 882 1155
Fax +1 516 467 3125
Web Site <http://www.heka.com>
Email ussales@heka.com
support@heka.com

General notice:

Product names used herein are for identification purposes only and may be trademarks of their respective owners. HEKA disclaims any and all rights in those marks.

We reserve the right to effect technical changes as development progresses. Special versions are available on request. Further technical data are provided by a detailed description, which is available on request. A guarantee of one year applies on all instruments.